

lifetime limited warranty

Your Drive brand product is warranted to be free of defects in materials and workmanship for the lifetime of the original consumer purchaser.

This device was built to exacting standards and carefully inspected prior to shipment. This Lifetime Limited Warranty is an expression of our confidence in the materials and workmanship of our products and our assurance to the consumer of years of dependable service.

In the event of a defect covered by this warranty, we will, at our option, repair or replace the device.

This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, casters, and grips, which are subject to normal wear and need periodic replacement.

If you have a question about your Drive device or this warranty, please contact an authorized Drive dealer.



© 2006 Medical Depot, Inc. All rights reserved.

Drive is a trademark of Medical Depot, Inc.

Port Washington N.Y. 11050 USA Made in China

super heavy duty bath bench

drive
MEDICAL DESIGN & MANUFACTURING

www.drivemedical.com

super heavy duty bath bench

Instructions for item # DR12029-2.



Item # DR12029-2

other bariatric products available



HeAvy Duty
TrAnsfer BenCh
Item # DR12025



Oversized Commode
Item # DR11117



Super Heavy Duty
Transfer Bench
Item # DR12009-2



Super Heavy Duty
Extra Large
Drop Arm Commode
Item # DR11132-1

caring for your bath bench

Never use bath safety seat if any looseness in the assembly is noticed, as personal injury may result. Sit, do not stand on the seat to avoid possible dangerous fall.

CAUTION: Maximum weight limitations for the product is not to exceed 600 lbs. Meets ANSI Standard A1171 Specifications.

NOTE: To maintain the beauty of the original finish, clean with mild soap and rinse with warm water. Wipe bath seat dry after each use.